

1. Product Information

| Product name | Recombinant spCas9 | |
|---------------------------------|--|--|
| Internal product code | BMC9 | |
| Catalog numbers (aliquot sizes) | 540-025-GMP (2.5 mg) 540-050-GMP (5 mg) | |
| Quality grade | current Good Manufacturing Practice (cGMP) | |
| Description | Recombinant wild type Cas9 nuclease from <i>Streptococcus</i> pyogenes with nuclear localization sequence (NLS-spCas9-NLS) | |
| Molecular weight | 162 kDa | |
| Concentration | 10 mg/mL solution in buffer (frozen) | |
| Storage buffer | 25 mM Tris, 300 mM NaCl, 0.1 mM EDTA, 50% glycerol, pH 7.4 buffer solution | |
| Storage temperature | -75 ± 15 °C | |
| Lot. Nr. | Specified on product label | |
| Use by date | Specified on product label | |
| Manufacturer | Biomay AG, Ada-Lovelace-Straße 2; 1220 Vienna, Austria; www.biomay.com; info@biomay.com | |

2. Description

CRISPR nuclease Cas9 (spCas9 from *Streptococcus pyogenes*, uniprot Q99ZW2 (CAS9_STRP1)) with nuclear localization sequences (NLS) on the N- and C-terminus. The product has been expressed as recombinant protein in *Escherichia coli*, purified by chromatography, filtered (0.2µm) and aseptically filled as a sterile product.

3. Intended Use / Application

cGMP grade: product is a recombinant protein that has been manufactured and quality-controlled under the conditions of current Good Manufacturing Practice (GMP). It has been certified and released by a Qualified Person (QP) under EMA law (EMA directive 2001/83/EC). It was designed and is intended to be used for gene-editing of eukaryotic cells with a specific guide RNA (gRNA). It is fully compliant to be used as critical starting material for gene-editing in clinical trials (phase I-III) or for commercial use.

Biomay imposes no restrictions on clients regarding the use of the purchased Cas9, except in the field of hemoglobinopathies. In the specific area of hemoglobinopathies, Cas9 obtained from Biomay must not be used for preclinical or clinical trials (Phases I, II or III), nor for any commercial applications. Clients are required to acknowledge and accept these specific terms related to Biomay's Cas9 with the acquisition of the material.

Version 2 (Sep 2025) Page 1 of 2



4. Quality Control and Specifications

| Parameter | Quality Control Method | Specification cGMP grade |
|--------------------------|--|--|
| Appearance | Visual inspection | Clear and colourless liquid, free of visible particles |
| pH-value | pH potentiometric (Ph. Eur. 2.2.3) | 7.4 ± 0.3 |
| Content (concentration) | UV spectrophotometric (UV 280 nm / Ph. Eur. 2.2.25) | 9.0 - 11.0 mg/ml |
| Homogeneity (aggregates) | Size exclusion chromatography (HP-SEC) | ≥ 95% monomer |
| Identity | Western blot | Positive reaction of the main band with the Cas9 specific antibody |
| | | Migration of the main band conforms to reference |
| Purity | SDS-PAGE | Migration of the main band conforms to reference |
| | | Purity: ≥ 95% |
| Host cell protein | E. <i>coli</i> host cell protein ELISA | ≤ 350 ng/mg (≤ 0.035 % w/w) |
| Endotoxins | LAL test (Ph. Eur. 2.6.14 method D) | < 5 EU/mg |
| Residual RNase | Fluorimetric | < 15 mU/ml |
| Residual DNase | Fluorimetric | No DNase activity detectable |
| Host cell DNA | qPCR | ≤ 30 ng/mg (≤ 0.003 % w/w) |
| Sterility | Sterility testing (Ph. Eur. 2.6.1) | Sterile |
| Functionality assay | RP-HPLC | 70 – 150% relative EC ₅₀ |

5. Safety information

Material is considered non-infectious, non-toxic and non-pathogenic under the conditions of the intended applications. General safety procedures should still be followed to maintain a safe working environment. Always wear appropriate personal protective equipment (PPE), including lab coats, gloves, and safety glasses, to avoid contamination and accidental exposure. Work in a clean, organized space, and handle reagents with care, avoiding direct contact. Dispose of all waste materials, including gloves and pipette tips, in designated biohazard containers, even though the commodity is non-hazardous, to prevent cross-contamination. Always wash hands after handling any biological material and before leaving the lab.

Version 2 (Sep 2025) Page 2 of 2