# Alt a 6 (formerly Alt a 2; Alt a 6.0101 according to effective allergen nomenclature)

(enolase from Alternaria alternata)

For research purpose only.

## PRODUCT DESCRIPTION:

Lot#: 02

Amount: 1 mg

Quality: Purity approx. 98%.

Access: EMBL: U82437 Uniprot: Q9HDT3

Mw = 47,382 Dalton

Mol. Ext. Coeff.: 48,485; 1 mg/ml A<sub>280</sub>=1.023\*

pl = 5.3

#### General information:

BIOMAY Alt a 2 (Alt a 6 according to the effective allergen nomenclature) is a lyophilized, recombinant protein with IgE-binding capacity. It was produced by heterologous expression in *E. coli*, purified by conventional biochemical methods and lyophilized from 5 mM sodium phosphate buffer pH 7.4.

#### Reconstitution:

The material can be reconstituted with distilled water (or equivalent) or dilute buffers. Addition of 1-2 mM  $\beta$ -Mercaptoethanol is recommended. If reconstituted to a concentration of 1 mg/ml, the sodium phosphate buffer concentration will be 2 mM (pH 7,4). The product is soluble to more than 95% at a concentration of 2 mg/ml.

### Storage:

When stored at -20°C the quality of the material will be maintained for several years. However, for short periods (max. 3 weeks) the lyophilized product may be kept at room temperature. After reconstitution store at -20°C. Avoid repeated freezing/thawing.

biomay

Vienna Competence Center Lazarettgasse 19 Top 1

Tel: +43 1 7966296-0

Fax: +43 1 7966296-111 e-mail: info@biomay.com www.biomay.com

BIOMAY AG

A-1090 Wien

#### Quality control:

Purity has been determined by SDS-PAGE and staining with Coomassie Brilliant Blue R-250.

Alt a 2 (Alt a 6) Lot# 02 tested positive in an IgE-Immunoblot with a standardized pool of human Alt a 2 (Alt a 6)-reactive sera.

Enolase activity can be restored by dissolving the lyophilized protein to a concentration of 1 mg/ml in 4.5 M urea and dialysing it over night against 10mM Na<sub>x</sub>H<sub>x</sub>PO<sub>4</sub> pH 7.5, containing 2mM  $\beta$ -Mercaptoethanole, and 2mM Mg<sup>2+</sup>.

No enzymatic activity was observed after reconstitution in water/ $\beta$ -Mercaptoethanole.

The mol.ext.coeff. was calculated from the DNA-derived protein sequence as described by Gill, S.C. and by Hippel, P.H. (1989), Analytical Biochemistry 182, 319-326.