# PRODUCT DESCRIPTION

# Der p 5 (Der p 5.0101)

(Dermatophagoides pteronyssinus, house dust mite allergen 5)

# For research purpose only.

Access: Genbank: X17699 UniProt: P14004 Molecular Weight: 13,718 Dalton (according to sequence from database without signal peptide but

with N-terminal methionine)

Mol. Ext. Coeff.: 4,470; 1 mg/ml A<sub>280</sub>=0.326

pl = 5.16



## **BIOMAY AG**

Vienna Competence Center Lazarettgasse 19 Top 1 A-1090 Wien

Tel: +43 1 7966296-0 Fax: +43 1 7966296-111

e-mail: info@biomay.com

Purity: ≥ 99% (SDS-PAGE) Endotoxin content: < 0.0005 EU/µg

IgE-Reactivity: reacts with IgE from Der p 5

reactive human serum

Lot#: 01

Amount: 1 mg

#### General information:

BIOMAY Der p 5.0101 is a recombinant protein with IgE-binding capacity. It was produced by heterologous expression in E. coli, purified by conventional biochemical methods, and lyophilized from 5 mM sodium phosphate buffer pH 7.4.

## Quality control of the product:

Purity has of the product been determined on SDS-PAGE gels stained with Coomassie Brilliant Blue R-250. Endotoxin content was determined by using a Limulus Amebocyte Lysate (LAL) assay. The lot stated above tested positive in an IgE-Immunoblot with Der p 5 -reactive serum.

# Storage of lyophilized product:

When stored at ≤ -20°C the quality of the lyophilized material is maintained for several years (see expiration date on the vial). For short periods (max. 3 weeks) the lyophilized product may be kept at room temperature.

## Reconstitution properties:

To achieve a solubilization of the product of at least 90%, we recommend to reconstitute the lyophilized protein to a concentration of 1 mg/mL with water of appropriate quality. After complete reconstitution the product concentration can be adjusted with the desired buffer as required, whereby the product must be principally soluble under the conditions applied.

# Reconstitution procedure:

Carefully inspect the vial for the location of the lyophilisate pellet. Some lyophilisates or pieces thereof are loose and might be located near the cap. In this case spin down the lyophlisate in a suitable centrifuge. Open the cap just as wide as necessary and pipet 1000 µL of water of appropriate quality into the vial. Close the cap and invert the vial several times, so that the complete lyophilisate and the whole inner surface of the vial are wetted. Incubate the vial for 2 h at room temperature on a rolling or an overhead incubator. Alternatively manual agitation can be applied by inverting the vial several times followed by gentle vortexing. This manual agitation procedure should be repeated several times during the incubation time. After the incubation time carefully visually inspect the tube for remaining undissolved material and eventually continue the incubation until the product is completely dissolved.

# Storage of reconstituted product:

Reconstituted product which is not used directly after reconstitution should be aliquoted in small aliquots (10-50 µL) and stored at ≤ -20°C. After thawing use these aliquots at once and avoid repeated freezing/thawing cycles.

The mol.ext. coeff. was calculated from the DNA-derived protein sequence as described by Gill, S.C. and by Hippel, P.H. (1989), Analytical Biochemistry 182, 319-326.

# PRODUCT DESCRIPTION

# Der p 5 (Der p 5.0101)

(Dermatophagoides pteronyssinus, house dust mite allergen 5)

# For research purpose only.

Access: Genbank: X17699 UniProt: P14004

Molecular Weight: 13,718 Dalton (according to sequence from database without signal peptide but

with N-terminal methionine)

Mol. Ext. Coeff.: 4,470; 1 mg/ml A<sub>280</sub>=0.326

pl = 5.16



## **BIOMAY AG**

Vienna Competence Center Lazarettgasse 19 Top 1 A-1090 Wien

Tel: +43 1 7966296-0 Fax: +43 1 7966296-111 e-mail: info@biomay.com

e-mail: info@biomay.c

Endotoxin content: < 0.0005 EU/µg

IgE-Reactivity: reacts with IgE from Der p 5

reactive human serum

Purity: ≥ 99% (SDS-PAGE)

Lot#: 01

Amount: 250 µg

## General information:

BIOMAY Der p 5.0101 is a recombinant protein with IgE-binding capacity. It was produced by heterologous expression in *E. coli*, purified by conventional biochemical methods, and lyophilized from 5 mM sodium phosphate buffer pH 7.4.

## Quality control of the product:

Purity has of the product been determined on SDS-PAGE gels stained with Coomassie Brilliant Blue R-250. Endotoxin content was determined by using a Limulus Amebocyte Lysate (LAL) assay. The lot stated above tested positive in an IgE-Immunoblot with Der p 5-reactive serum.

# Storage of lyophilized product:

When stored at  $\leq -20^{\circ}$ C the quality of the lyophilized material is maintained for several years (see expiration date on the vial). For short periods (max. 3 weeks) the lyophilized product may be kept at room temperature.

# Reconstitution properties:

To achieve a solubilization of the product of at least 90%, we recommend to reconstitute the lyophilized protein to a concentration of 1 mg/mL with water of appropriate quality. After complete reconstitution the product concentration can be adjusted with the desired buffer as required, whereby the product must be principally soluble under the conditions applied.

# Reconstitution procedure:

Carefully inspect the vial for the location of the lyophilisate pellet. Some lyophilisates or pieces thereof are loose and might be located near the cap. In this case spin down the lyophlisate in a suitable centrifuge. Open the cap just as wide as necessary and pipet 250 µL of water of appropriate quality into the vial. Close the cap and invert the vial several times, so that the complete lyophilisate and the whole inner surface of the vial are wetted. Incubate the vial for 2 h at room temperature on a rolling or an overhead incubator. Alternatively manual agitation can be applied by inverting the vial several times followed by gentle vortexing. This manual agitation procedure should be repeated several times during the incubation time. After the incubation time carefully visually inspect the tube for remaining undissolved material and eventually continue the incubation until the product is completely dissolved.

# Storage of reconstituted product:

Reconstituted product which is not used directly after reconstitution should be aliquoted in small aliquots (10-50  $\mu$ L) and stored at  $\leq$  -20°C. After thawing use these aliquots at once and avoid repeated freezing/thawing cycles.

The mol.ext. coeff. was calculated from the DNA-derived protein sequence as described by Gill, S.C. and by Hippel, P.H. (1989), Analytical Biochemistry <u>182</u>, 319-326.